PRODUCT INFORMATION

SERVA DNA Stain G

Cat. No.: 39803

PRODUCT DESCRIPTION

SERVA DNA Stain G is a nucleic acid stain which can be used as a safer alternative to the traditional ethidium bromide stain for detecting nucleic acid in agarose gels. It is as sensitive as ethidium bromide and can be used exactly the same way in agarose gel electrophoresis. SERVA DNA Stain G emits green fluorescence when bound to DNA or RNA. It has two fluorescence excitation maxima when bound to nucleic acid:

One at approx. 300 nm and another at approx. 450 nm.

1 ml of SERVA DNA Stain G is sufficient for 20 L of agarose gel.

Protocol:

- Prepare 100 ml of agarose gel solution (concentration from 0.8 - 3.0 %) and heat until the solution is completely clear and no small floating particles are visible.
- Add 2 - 5 µl of SERVA DNA Stain G to the gel solution and mix it gently.
- Cool the gel to 50 – 60 °C and cast the gel into the gel tray.
- When the gel is solid, load the samples and perform electrophoresis.
- Detect the bands under UV illuminator.

Notes:

- The thickness of gel should < 0.5cm.
 - Repeated melting of gels containing SERVA DNA Stain G may result in low sensitivity.
 - SERVA DNA Stain G allows visualization of DNA (> 50 ng) in the agarose gel under visible light. This eliminates the need for exposure to UV light, which can nick and damage DNA.
 - SERVA DNA Stain G is non-carcinogenic but may irritate skin and eyes. Please wear gloves while handling.
-

Storage: Store at room temperature or at + 4 °C, protected from light.

Applications: Non-carcinogenic alternative to Ethidium bromide.

Safety: SERVA DNA Stain G is non-carcinogenic and according to the Ames test it causes significantly fewer mutations than ethidium bromide.

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